

GAMMA IRRADIATED

Klercide™ WFI Quality Water

For dilution and residue removal

Maintaining high quality throughout the cleanroom

Klercide WFI Quality Water is suitable for use in product manufacturing areas. It is produced from USP Water for Injection (WFI) with guaranteed endotoxin levels below 0.25 EU/mL.

Klercide WFI Quality Water is useful for the dilution of sterile concentrates to maintain the required quality of diluted product for high grade cleanrooms.

A highly effective way of removing any residues after disinfection, simply wipe down with Klercide WFI Quality Water after using the Klercide range of biocides in areas where residues may be an issue.

Klercide WFI Quality Water is also:

- Available double or triple bagged
- Supplied in our patented and validated SteriShield™ Delivery System (SDS)

The SteriShield Delivery System (SDS) features:

- A closed trigger system providing a completely sealed system
- Protection for the integrity of the contents throughout use
- Validated in-use shelf life of three months
- A fully adjustable spray pattern to ensure optimum surface coverage



ECOLAB®

GLOBAL EXPERTS IN
CONTAMINATION CONTROL

Why use Klercide WFI Quality Water?

Effective

Klercide WFI Quality Water is particularly effective at removing undesirable compounds and diluting concentrated products.

It is 0.2 micron filtered, filled and capped in an ISO Class 7 cleanroom.

Cost effective

The product is available in either a 32 oz trigger spray or a 1 US gallon capped container. The 32 oz trigger spray is fitted with the SteriShield Deliver System (SDS), which ensures all contents are dispensed from the bottle so there is no wastage. The SDS has been validated to enable an in-use shelf life of three months so there is no need to discard unused product.

Designed for cleanrooms

Klercide WFI Quality Water is manufactured in Ecolab's ISO 9001 accredited facilities to the principles of cGMP. These manufacturing conditions ensure Klercide WFI Quality Water is manufactured to the highest quality and are compliant for use in your critical areas.

Klercide WFI Quality Water is:

- ▶ Endotoxin tested to an approved pharmacopoeia method
- ▶ Supported with Certificates of Analysis and Sterility with every batch
- ▶ Supported with full documentation



Klercide WFI Quality Water
32 oz SDS spray



Klercide WFI Quality Water
1 US gal capped

For more information contact your Ecolab Account Manager

WORLDWIDE HEADQUARTERS

1 Ecolab Place
St. Paul, MN 55102
USA

www.ecolab.com/lifesciences

GAMMA IRRADIATED

Klercide™ WFI Quality Water

For dilution and residue removal

Technical specifications

Application area

For the preparation or dilution of disinfectants or for residue removal.

Instructions for use

For use instructions, please see the product label. Not for parenteral administration.

Packaging information

Klercide WFI Quality Water

Unit size	Unit of sale
32 oz spray	6 x 32 oz
1 US gal capped	4 x 1 US gal

For product order codes and additional information, please contact your account manager or visit our website

Shelf life

Unopened, this product has a two year shelf life.

USE BIOCIDES SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

