

Ecolab[®] Disinfectant 1 Wipe With Plastic-Free Wipes¹

Training Resource Summary

Preparing for Trial and Implementation

A What to know about the wipe experience.

- Wipes are made from **wood-pulp fibers**, making them **highly absorbent** with a **cloth-like feel**.
- **Fold larger wipes** before use for great handling.
- After removing a wipe from the roll, **small fibers at the end of the wipe may be noticed. This is normal** —use as usual and avoid intentionally pulling on the fibers.²

B Access training materials.

- **Video Trainings:** click the link or visit YouTube for Refillable Bucket & Canister Training In-service videos.
- **Guides and claims:** Refillable bucket training guide, end-user wipe tip sheet, kill claims summary, etc., see Ecolab website.

Click or visit:
[Ecolab.com/PFD1Wipe](https://www.ecolab.com/PFD1Wipe)



1. Wipe material is plastic-free (packaging made of plastic: bucket, refill bag). EPA Reg. No. 1677-263. Kills 40+ organisms in 60 seconds and kills *C. auris* and *T. interdigitale* in 75 seconds. Refer to product label for full organism listing & use directions.
2. See FAQ on website for more information.