

Sink & Surface Cleaner Sanitizer Test Strips How-To Guide

WHAT YOU WILL NEED:

Sink & Surface Cleaner Sanitizer Test Strip Vial

Item Number: 20318-01-11



HOW TO USE THE TEST STRIPS:

1. Withdraw a test strip from the canister.



2. Dip Test Strip for 5 seconds into the test solution*. Shake off excess solution.



3. Compare colors after 10 seconds with colors on the test strip canister to determine the concentration.



* The tests strips will only work with Sink & Surface Cleaner Sanitizer and cannot be used with other products. The product has 2 active ingredients; the test strips are designed to read both active ingredients at the same time.

HOW TO READ THE LABEL:

The first row of numbers indicates the use dilution of Sink & Surface Cleaner Sanitizer in oz/gal.

The approved dilution range of the product is 0.27 - 0.55 oz/gal.

The blue box indicates the EPA approved active range of Sink & Surface Cleaner Sanitizer. The color of the test strip should fall within the color range indicated here.

The second row of numbers indicates the use dilution of Sink & Surface Cleaner Sanitizer in mL/L.

The approved dilution range of the product is 2.11 - 4.30 mL/L.

The third row of numbers indicates the level of DDBSA (one of the two approved actives - the test strips are designed to read both active ingredients at the same time) in the use dilution of Sink & Surface Cleaner Sanitizer.

The approved active range of DDBSA is 272 - 700 ppm.

both active ingredients at the same time) in the use dilution of Sink & Surface Cleaner Sanitizer.

The approved active range of Lactic Acid is 704 - 1875 ppm.

ounces / gallon		ounces / gallon		onzas / galón	
0	0.15	0.27	0.55	0.75	1.00
milliliters / liter		millilitres / litre		mililitros / litro	
0	1.17	2.11	4.30	5.86	7.81
					
0	170	272	700	848	1130
ppm DDBSA*					
0	452	704	1875	2258	3010
ppm Lactic Acid					

Contact your Ecolab Representative or call 1 800 352 5326 for more information.